

The StairWalker

Blue A

Walking Aids with Stairs

UCSF Disabilities Statistics Center study on mobility device use:

- **1.8 million** Americans use walkers
- Only 39.3% of mobility device users live in houses with just one floor
- 62.1% of mobility device users must climb steps to get in or out of their homes
- Only 9.6% of mobility device users have added stair-climbing features

Our contact, Luis Loya (Access Planner for Kessler McGuinness and Associates) verifies need for product.



Current Products

Stair lifts:

- Person sits on chair, which then climbs up the side of the wall



Expensive:

- Cheapest are about \$2,500


No walker:

- Will not bring up the walker, so person lands on the upper floor without walking aid

Current Products

US Patent 6,453,921B1:

- Retractable legs
- Activated by gripping handle




Safety risk → **falling**

- 1/3 of people 65 and older fall at least once a year
- Fear of falling is main fear among elderly, above robbery and financial fears


Requires **hand dexterity**, which elderly often don't have

The StairWalker

- Attaches to the railings
- As the user pushes the walker up against the bar, the StairWalker moves up too
- The StairWalker will not fall back down



How It Works



- Walker attaches to a crossbar at the front
 - Rubber rollers to firmly grip the railing
 - Uses **ratchets** to keep from falling down

Features

- **Safety** – prevents user from falling
- **Mobility** – helps users have full access to their walkers
- Low **cost**
- Requires minimal **installation**


Critical parts:

- Properly compressible rubber rollers to grip rail
- Sliding enough to push it up
- Interfacing with the wall and existing railing

Further Considerations

Range of railing shapes and sizes:

- Attach a baseline track to the railing, sold with the StairWalker.
- The box could clamp onto the railing from the top, or from the side →



Going down stairs:


- Orient the rollers on each side in different directions. Lock one side each time, depending on direction.

For railings that do not start at the first step, include rail extension.

Further Considerations

Range of railing shapes and sizes:

- Attach a baseline track to the railing, sold with the StairWalker.
- The box could clamp onto the railing from the top, or from the side →



Going down stairs:

- Orient the rollers on each side in different directions. Lock one side each time, depending on direction.

For railings that do not start at the first step, include rail extension.

StairWalker

Can be an aid and asset to the 1.8 million Americans who use walkers.

References

- Kaye, H. Stephen, Taewoon Kang, Mitchell P. LaPlante. “Mobility Device Use in the United States.” *UCSF Disabilities Statistics Center*, June 2000
http://dsc.ucsf.edu/publication.php?pub_id=2§ion_id=1 (1 October 2008).
- “Use walkers to lower the risk of falling” *India eNews*, 5 February 2008
<http://www.indiaenews.com/europe/20080205/95666.htm> (1 October 2008).
- Howland, Jonathan, et al. “Fear of Falling among the Community-Dwelling Elderly”, *Journal of Aging and Health*, 5, No. 2 (1993): 229-243.
- Rost, Brian M. “Stair walker.” US Patent 6,453,921 B1. 24 September 2002.